


PRODUCTS CATALOG



Quality is our top priority

LABORATOIRES
ellios
B i o T e k

25 of innovation in design
years of Disposable Medical Devices
for Assisted Reproductive Technology

TABLE OF CONTENTS

ASSISTED REPRODUCTION

• ELLIOCATH N.V.[®] – Embryo transfer catheter	p. 5
• ELLIOCATH S[®] – Echogenic embryo transfer catheter	p. 7
• ELLIOSET – Echogenic embryo transfer set	p. 9
• EMBRYO TRANSFER SET – With removable mandrel – For difficult transfer	p. 11
• OOCYTES ASPIRATION NEEDLE – Back cut bevel	p. 13
• OOCYTES PUNCTURE SET – For Falcon 14 mL	p. 15
• OOCYTES PUNCTURE SET – For Falcon 50 mL	p. 17
• OOCYTES PUNCTURE SET – For vial sv-XL30	p. 19
• VIAL SV-XL30 – For oocytes puncture set with screw cap	p. 21
• VACUUM LINE WITH HYDROPHOBIC – Filter	p. 23

FUNCTIONAL EXPLORATION

• HYSOKAT – For hysterasonography - Hystero-salpingo-foam-sonographie (HyFoSy)	p. 25
• SONOKAT – For hysterasonography/insemination	p. 27



Reference	Product	Sales Unit	I.D.	Length	Sterilisation
EBT018-KT-NV	ELLIOCATH N.V.®	25 units Individually packed	1.1 mm	20.5 cm	γ irradiation

ELLIOCATH N.V.®

Embryo transfer catheter (metal free)

PRODUCT DESCRIPTION

- Curved transfer catheter making easier the passage through the cervical canal.
- Tactile marks on the grip of ELLIOCATH N.V.® indicate the direction of the curve of the catheter.
- Set of 4 labels for traceability.
- Centimetric graduation on 360°.
- Tight sliding ring.

INDICATIONS

- Embryo transfer.
- Intra Uterine Insemination.

RECOMMENDATIONS FOR USE


- Before use, wash carefully the catheter with an appropriate sterile culture medium.

Embryo transfer

- Place the embryo in ELLIOCATH N.V.® with a small volume of culture medium.
- Position the sliding ring according to the uterus depth previously measured.
- Advance ELLIOCATH N.V.® to a distance 1 cm from the uterine fundus.
- To facilitate cervical canal passage, push the catheter with a gentle rotational movement.
- Transfer embryo.
- Slowly withdraw the catheter.

Intra uterine insemination

- Fill ELLIOCATH N.V.®, connected to a syringe, with the sperm preparation.
- Push ELLIOCATH N.V.® into the uterus.
- To facilitate cervical canal passage, push the catheter with a gentle rotational movement.
- Push slowly the sperm preparation.
- Slowly withdraw the catheter.



Reference	Product	Sales Unit	I.D.	Length	Sterilisation
EBT018-KT-S	ELLIOCATH S®	25 units Individually packed	1.1 mm	20.5 cm	γ irradiation

ELLIOCATH S®

Echogenic embryo transfer catheter

PRODUCT DESCRIPTION


- Precurved transfer catheter, with an atraumatic echogenic bulb tip:
 - making easier passage through cervix,
 - providing easy identification of its position.
- Tactile marks on the grip indicating the direction of the curve of the catheter.
- Adjustable positioner atraumatic easy to handle.
- Set of 4 labels for traceability.

INDICATIONS

- Embryo transfer.

RECOMMENDATIONS FOR USE

- Before use, wash carefully the catheter with an appropriate sterile culture medium.
- Place the embryo in ELLIOCATH S® with a small volume of culture medium.
- Position the sliding ring according to the uterus depth previously measured.
- Advance ELLIOCATH S® to a distance 1 cm from the uterine fundus.
- To facilitate cervical canal passage, push the catheter with a gentle rotational movement.
- Transfer embryo.
- Slowly withdraw the catheter.



References	Product	Sales Unit	Specifications			Sterilisation
EBT023-ST-S	Echogenic embryo transfer set	10 units Individually packed		Guiding catheter	Transfer Catheter	γ irradiation
			I.D.	1.7	0.7	
			O.D.	2.2	1.5	
			Transfer catheter exceeds over 5 cm of the guiding catheter			

ELLIOSET

Echogenic embryo transfer set

PRODUCT DESCRIPTION

Guide catheter (soft P.P.)

- Precurved, with an atraumatic **echogenic bulb tip**:
 - making easier passage through cervix,
 - providing easy identification of its position and atraumatic entry to cervical canal.
- Adjustable positioner atraumatic easy to handle.
- Grip shape facilitates introduction of the transfer catheter.

Transfer catheter (soft P.U.)

- **Echogenic tip.**

Centimetric marks (on guide and transfer catheters) for confident placement of embryos.


Set of 4 labels for traceability.

INDICATIONS

- Embryo transfer.

RECOMMENDATIONS FOR USE

- Before use, rinse carefully the embryo transfer catheter with an appropriate sterile culture medium.
- Once embryos transferred, withdraw simultaneously and very slowly the 2 catheters.



References	Product	Sales unit	Specifications			Sterilisation
EBT00118-3P	Embryo transfer set 18 cm	10 units Individually packed		Guiding catheter	Transfer Catheter	γ irradiation
			I.D.	1.7	0.7	
			O.D.	2.3	1.5	
			Transfer catheter exceeds over 5 cm of the guiding catheter			
EBT00123-3P	Embryo transfer set 23 cm	10 units Individually packed		Guiding catheter	Transfer Catheter	γ irradiation
			I.D.	1.7	0.7	
			O.D.	2.3	1.5	
			Transfer catheter exceeds over 5 cm of the guiding catheter			

EMBRYO TRANSFER SET

Set for difficult embryo transfer

PRODUCT DESCRIPTION


- Guided catheter (polypropylene), with a sliding ring with centimetric graduations.
- Malleable removable mandrel.
- Transfer catheter (polyurethane) with centimetric graduations on the proximal part.
- **The proximal (luer) end of the guiding catheter has been specifically designed for making easier the transfer catheter insertion.**

INDICATIONS

- Embryo transfer in uterus, especially in the case of difficult transfers.
- Intra Uterine Insemination.

RECOMMENDATIONS FOR USE

- Before use, wash carefully the embryo transfer set with an appropriate sterile culture medium.
- With the help of centimetric graduations, position the sliding ring (guiding catheter) according to the uterus depth previously measured.
- Give a shape if necessary to the guiding catheter + mandrel. Cross the cervical canal up until the internal orifice.
- Withdraw the mandrel.
- Insert the transfer catheter into the guiding catheter and push it up until the first centimetric graduation reaches the guiding catheter entrance. Distal ends of both catheters match properly.
- Then push the transfer catheter with the help of centimetric graduations and proceed to the embryo transfer.
- Withdraw simultaneously and very slowly the 2 catheters.



Reference	Product	Sales unit	O.D.	Length	Sterilisation
EBT00030-BC	Ovum aspiration needle with stylet	10 units individually packed	17 G (1.4 mm)	30 cm (active length)	γ irradiation
EBT00030-BC-0	Ovum aspiration needle without stylet	10 units individually packed	17 G (1.4 mm)	30 cm (active length)	γ irradiation

OOCYTES ASPIRATION NEEDLE

Back cut bevel 17 G (=1.4 mm)

PRODUCT DESCRIPTION

- Single lumen needle.
- Back-cut bevel.
- Non-Coring bevel.
- 2 laser echo markings.
- Plastic stylet.
- Large ergonomic handle.
- Tactile mark indicating bevel position.
- 4 tabs for traceability.

INDICATIONS

- Echoguided oocytes pick up.
- Ovarian cyst aspiration.

RECOMMANDATIONS FOR USE

- Before use, wash carefully the needle with an appropriate sterile culture medium.
- Adapt a syringe to the luer tip.
- Proceed according to the selected route.



References	Products	Sales unit	Specifications	Sterilisation
EBT01230S-0	Follicule aspiration set (Falcon 14 mL)	10 units individually packed	Needle: 30 cm - 17 G (1.4 mm) Needle – stopper : 90 cm Stopper – pump: 100 cm	Ethylene oxide
EBT00230S-0	Follicule aspiration set (Falcon 14 mL)	10 units individually packed	Needle: 30 cm - 17 G (1.4 mm) Needle – stopper : 90 cm Stopper – pump: 100 cm	γ irradiation
EBT00230S-R	Follicule aspiration set with stopcock (Falcon 14 mL)	10 units individually packed	Needle: 30 cm - 17 G (1.4 mm) Needle – stopper : 90 cm Stopper – pump: 100 cm	γ irradiation
EBT00235S-0	Follicule aspiration set (Falcon 14 mL)	10 units individually packed	Needle: 35 cm - 17 G (1.4 mm) Needle – stopper: 80 cm Stopper – pump: 90 cm	γ irradiation
EBT00235S-R	Follicule aspiration set with stopcock (Falcon 14 mL)	10 units individually packed	Needle: 35 cm - 17 G (1.4 mm) Needle – stopper : 80 cm Stopper – pump: 90 cm	γ irradiation

OOCYTES PUNCTURE SET

for Falcon 14 mL

Back cut bevel 17 G (=1.4 mm)

PRODUCT DESCRIPTION

- Single lumen needle.
- Back-cut bevel.
- Non-Coring bevel (prevents tissue entering lumen).
- 2 laser echo markings.
- Tactile mark indicating bevel position.
- 2 tubes connected to a silicone stopper for 14 mL Falcon-like tube.
- 4 tabs for traceability.


INDICATIONS

- Echoguided oocytes pick up.
- Ovarian Cyst aspiration.

RECOMMANDATIONS FOR USE

- Before use, wash carefully the needle and the aspiration line with a appropriate sterile culture medium.
- Connect a first Falcon-type sterile tube to the stopper (to be replaced when needed during the procedure).
- Proceed to the ovum aspiration under controlled pressure.
- Rinse out carefully once the aspiration is ended (oocytes, if present, will be collected from the needle and the aspiration line).

The oocytes puncture set with a 3 ways stopcock allows intermittent flushing and aspiration.



References	Products	Sales unit	Specifications	Sterilisation
EBT00430S-0	Follicule aspiration set (Falcon 50 mL)	10 units individually packed	Needle: 30 cm - 17 G (1.4 mm) Needle – stopper : 90 cm Stopper – pump: 100 cm	γ irradiation

SET DE PONCTION D'OVOCYTES

pour Falcon 50 mL

Aiguille back cut 17 G (diam. 1.4 mm)

PRODUCT DESCRIPTION


- Single lumen needle.
- Back-cut bevel.
- Non-Coring bevel (prevents tissue entering lumen).
- 2 laser echo markings.
- Tactile mark indicating bevel position.
- 2 tubes connected to a silicone stopper for 50 mL Falcon-like tube.
- 4 tabs for traceability.

INDICATIONS

- Echoguided oocytes pick up.
- Ovarian Cyst aspiration.

RECOMMENDATIONS FOR USE

- Before use, wash carefully the needle and the aspiration line with a appropriate sterile culture medium.
- Connect a first Falcon-type sterile tube to the stopper (to be replaced when needed during the procedure).
- Proceed to the ovum aspiration under controlled pressure.
- Rinse out carefully once the aspiration is ended (oocytes, if present, will be collected from the needle and the aspiration line).



References	Products	Sales unit	Specifications	Sterilisation
EBT00330S-0	Follicule aspiration set with screw cap	10 units individually packed	Needle: 30 cm - 17 G (1.4 mm) Needle – cap: 90 cm Cap – pump: 100 cm	γ irradiation
EBT00330S-R	Follicule aspiration set with screw cap and with stopcock	10 units individually packed	Needle: 30 cm - 17 G (1.4 mm) Needle – cap: 90 cm Cap – pump: 100 cm	γ irradiation

OOCYTES PUNCTURE SET

for vial sv-XL30

Back cut bevel 17 G (=1.4 mm)

PRODUCT DESCRIPTION

- Single lumen needle.
- Back-cut bevel.
- Non-Coring bevel (prevents tissue entering lumen).
- 2 laser echo markings.
- Tactile mark indicating bevel position.
- 2 tubes connected to a screw cap (to be used with a 30 mL vial - sv-XL30).
- 4 tabs for traceability.

INDICATION

- Echoguided oocytes pick up.

RECOMMANDATIONS FOR USE

- Before use, wash carefully the needle and the aspiration line with an appropriate sterile culture medium.
- Screw to a vial sv-XL30 (to be replaced when needed during the procedure).
- Proceed to the ovum aspiration under controlled pressure.
- Rinse out carefully once the aspiration is ended (oocytes, if present, will be collected from the needle and the aspiration line).

The oocytes puncture set with a 3 ways stopcock allows intermittent flushing and aspiration.



Reference	Product	Sales unit	Sterilisation
EBT003050BV	Vial sv-XL30 for oocytes collection	50 units individually packed	γ irradiation

VIAL sv-XL30

for oocytes puncture set with screw cap

PRODUCT DESCRIPTION

- Specially designed to maintain follicular fluid temperature. When the vial is placed within a heat-controlled nest (sold by Carteau & Ost, Labotect, etc...) the flange located at the upper part (diameter equal to the vial base) prevents outside air arrival potentially leading to turmoil and heat shocks detrimental to oocytes.
- Cap and vial components selected to avoid any deformation at 37°C.
- Cork designed to prevent leakings.
- Handling easier with a quarter-turn screwthread.



References	Products	Sales unit	Specifications	Sterilisation
EBT01615	Hydrophobic filter for puncture set	25 units individually packed	Tube: 150 cm	Ethylene oxide

VACUUM LINE WITH HYDROPHOBIC Filter for puncture set

PRODUCT DESCRIPTION


- Hydrophobic filter 0.2 μm .
- Tubing of 150 cm with male ends connector.

INDICATIONS

- Oocytes puncture with aspiration pump.
- Match with all aspiration oocytes pump.

Biological risk

Never perform a follicular puncture if the follicular suction pump, tubing or filter shows signs of contamination. The tubing and its filter are consumables for single use. They must not be reused or resterilized. Reuse can cause cross-contamination leading to the transmission of infectious germs. Re-sterilization can damage the structure and performance of equipment that impacts the integrity of biological samples or suction equipment.



References	Product	Sales Unit	Length	Sterilisation
EBT01700	HYSOKAT	25 units individually packed	24.5 cm	Ethylene oxide

HYSOKAT

Hysterosonography - Hystero-salpingo-foam-sonographie (HyFoSy)

PRODUCT DESCRIPTION


- HYSOKAT features specially developed to perform Hystero-salpingo-foam sonography or traditional hysterosonography:
 - Central pre-curved soft catheter displaying:
 - Centimetric graduations,
 - Tactile marks on the grip indicating the direction of the curve.
 - External sliding tube providing:
 - A **very soft and atraumatic** cone at the distal end, for an effective and harmless occlusion of the cervix,
 - A jointly fixed ring, at the proximal end for precise positioning.

INDICATIONS

- Hysterosonography - Hystero-salpingo-foam-sonographie (HyFoSy).

RECOMMENDATIONS FOR USE

- Gently insert HYSOKAT through the cervix into the uterus.
- Push the external tube forward until an effective occlusion is obtained.
- Slowly inject the selected echogenic preparation.
- Upon completion of the procedure, gently withdraw HYSOKAT and properly discard according to clinical waste treatment procedures.



Reference	Product	Sales unit	Lenght (mm)	Sterilisation
EBT00500	SONOKAT	10 units individually packed	225	γ irradiation

SONOKAT

Hysterosonography

PRODUCT DESCRIPTION

- An internal PEBD catheter.
- An obturating disk.
- A stiff external tube.
- A two ways stopcock.

The obturator allows an effective and atraumatic occlusion of external cervical os whatever is the anatomic diameters.

This device is sterilized by gamma radiation.

INDICATIONS

- Hyterosonography.

RECOMMANDATIONS FOR USE

- Insert SONOKAT through the cervix into the uterus and push it till an effective occlusion is obtained.
- Maintaining the catheter in position, inject slowly the preparation (saline water).
- Shut off the stopcock for the time necessary for the related procedure.
- Remove slowly the device after completion of the procedure.



Laboratoires Ellios Bio Tek – 3 rue de l'arrivée – Tour C.I.T. – 75015 PARIS
Tél. : +33 (0)1 44 26 22 80 - www.ellios-biotek.com